FIGURE 1		
ESPAÑOL (ES)	DEUTSCH (DE)	FRANÇAIS (FR)
A ADAPTADOR MASCARILLA	MASKENADAPTER	ADAPTEUR DE MASQUE
B PIEZA FRONTAL	GESICHTSTEIL	PIÈCE FACIALE
C TOMA MUESTREO	PROBENABNAHMEANSCHLUSS	ORIFICE D'ÉCHANTILLONNAGE
D TOMA GIRATORIA	DREHGELENKANSCHLUSS	RACCORD ARTICULÉ
E RECTA	GERADER TYP	MODÈLE DROIT
F TOMA GIRATORIA	DREHGELENKANSCHLUSS	RACCORD ARTICULÉ
G ACODADA	GEBOGENER TYP	TYP COUDÉ
H CONEXIÓN COFIA	KOPFHAUBEN-BEFESTIGUNG	FENTES POUR ATTACHES DU HARNAIS

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US Pat No's D573,709 & 8,136,523 and
EC Req'd Design No. 000780556-001
ISO 9001 / ISO 13485

©0086 © 2015 Hans Rudolph, inc.

ENGLISH (en) - Product Description - Hans Rudolph 7500 Series V2 Mask - Oro-Nasal NIV Reusable - With NO Vents and NO AAV

Intended Use:
The 7500 V2 Mask is a reusable, multi-patient, multi-use, Oro-Nasal NIV mask intended to provide a patient interface for applications of noninvasive ventilation. The mask is intended for use on adults (>30 kg) in homes, hospitals and other clinical settings as an accessory to ventilators which have adequate alarms and safety systems for ventilator failure and which are intended to administer positive pressure ventilation for treatment of respiratory failure or respiratory insufficiency.

Caution:
Federal law restricts this device to sale by or on physician's order.

Warnings:

- Shipped clean, nonsterile. If disinfected or sterile use is required, follow the procedures in this document prior to use of the mask.
- To avoid skin irritation, never use any cleaning solutions or chemicals on this mask or headgear other than mild detergent and water.
- Patients with facial hair may experience mask leakage.
- Do not block mask or patient circuit flow.
- This mask does not include an Anti-Asphyxia Valve (AAV) that would allow patient to breathe if ventilator failed.
- If supplemental oxygen is used, inhaled percent oxygen will vary depending upon pressure settings, patient breathing pattern, mask size & leak rate.
- Oxygen flow must be turned off when the ventilation system is not operating.
- Combination of this mask with devices other than for which mask is intended can alter mask safety and performance becoming unsafe or life-threatening.

Contraindications:

- Open wounds prone to infection.
- Hemodynamic or cardiorespiratory instability.
- Unconsciousness or patients unable to remove mask.
- Claustrophobia, anxiety, or other discomfort with Oro-Nasal mask.
- Facial or nasopharyngeal deformity or other inability to fit mask & seal properly.
- Excessive reflux, GI bleed or other conditions.
- Impaired cough reflex, hiatal hernia, or inability to swallow or clear secretions.
- Barotrauma, upper airway obstruction or facial trauma.
- Recent facial, esophageal, or gastric surgery.
- Patients under medication with drug that may cause vomiting.
- Patients requiring immediate intubation.

Risks and Complications:

- Biocompatibility to device materials (irritation, sensitivity, or other reaction).
- Irritation/abrasion/ulceration/pain/deformity/irritation (face, head, neck, eyes, mucosa, dental).
- Inhaled, aspirated, applied or handled contaminants (particulates, microbes, body fluids, chemicals).
- Infection (skin, eyes, mucosa, respiratory).
- Excessive rebreathed dead space or breathing resistance (dyspnea, atelectasis, barotrauma).
- Excessive/damaging sound level.
- Low secretion clearance (especially during upper respiratory infections).
- Inadequate or excessive mask operating pressure (barotrauma, dyspnea, atelectasis, gastric distention, flatulence).
- Adverse interactions with attached non-HRI accessories.

Mask Components (Fig 1):

- Mask face piece has plastic braces with slot openings for attachment of headgear strap mounting quick-release clips.
- Swivel Port assembly (Elbow or Straight Type) comprises sampling port, 22 mm (ISO 5356-1) swivel port.

Application & Fitting (Fig 3):

- Determine patient mask size with optional Mask Sizing Gauge (Fig 2).
- Place mask over patient's nose and mouth (1). Chin should fit into mask chin cup portion with top of mask sealing area on the bridge of the nose.
- Slide headgear over patient's head (2) with one of the lower quick release strap clips disconnected. Reconnect headgear strap clip after headgear is placed (3). Bottom straps should be positioned below ears and top straps above ears & below eyes.
- Adjust headgear Velcro™ strap tension by pulling back slightly on both the top and bottom straps (4 & 5).
- Adjust straps at crown of head by pulling each in opposing directions until reaching required tension and then reconnecting Velcro™ (6).
- Connect mask to patient circuit by following ventilator operating instructions.
- Turn on ventilator. If you detect leaks reposition mask or adjust strap tension. If leak continues, regardless of adjustments, try another mask size.

Safety Features Verification:

- Mask Removal with Quick-Release Headgear (Fig 4). To remove mask slide one of the lower strap quick-release clips slightly forward towards mask front. It will disconnect from mask; pull mask and headgear to opposite side or up over patient's head.

Disassembly for Cleaning:

- Separate headgear from mask.
- Separate swivel port assembly from mask and rubber cap plug from mask adapter sampling port.

Cleaning Mask face piece, swivel port components & headgear:

Mask & swivel port assembly:

- Soak components for 5 minutes in warm water with mild detergent (neutral pH).
- Hand wash components with sponge or soft brush; then rinse with warm tap water.
- Place components in warm water bath and agitate for two minutes.
- Rinse in clean tap water for at least 1 minute.
- Allow components to air dry or dry with clean, lint free cloth.

Headgear:

- Submerge headgear in warm soapy tap water; gently rub all areas.
- Rinse in clean tap water for 1 minute or until all soap is removed and air dry.
- Precaution: Do not use bleach, chlorine or alcohol based solutions to clean mask or headgear. These solutions can damage components. Direct sunlight exposure on mask and headgear components can cause deterioration and reduce product life.

High Level Disinfection (Mask and Swivel port components less Headgear)

- Thoroughly clean mask prior to disinfection or sterilization.
- For disinfection use CideX™ liquid glutaraldehyde or hot water pasteurization. **Liquid Chemical Disinfection:** Soak the device components in the liquid chemical solution per the manufacturer's instructions. Remove the components and soak in 1500 ml of sterile water for a minimum of 1 minute. Repeat step 2 and then dry the device with a clean lint free cloth (preferably sterile). **Hot Water Pasteurization:** Immerse all device component surfaces in hot water temperatures 71.1°C – 76.6°C for 30 minutes. Dry the device with a clean lint free cloth (preferably sterile).

Visually inspect all components for cleanliness, function and defects. If there are any residues, stains or organic debris then repeat the previous steps. Dispose of and replace all defective parts.

Reassemble:

- Press rubber cap plug into sampling port of mask adapter.
- Install swivel port assembly into face mask groove in opening. The flange of the mask adapter mates with the groove of the mask opening.
- Mount headgear to face piece with the snap in quick-release clips.

Functional Check:

- Swivel port joints swivel freely; swivel port assembly is completely engaged in mask face piece. Do not use if it is non-functional.

Storage

- Components should be completely dry prior to storage in a sealed poly bag.
- Label poly bag with disinfection/clean status, mask description, date and initials.

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